Part No.

Description

Use:

Î Ï Í Á5020T01

Pump Sampler Set (PTFE)

Reusable liquid pump





Each Set Comprises of a Carrying Case containing:

Part No.	Description	Quantity	Material
675 5020A09	Hand Pump	1	Plastic with polymer O-ring seals
675 5020A11	Bottle Adapter*	1	PTFE
675 5020T05	Sampling Tube (OD 8mm) x 5m length)	1	PTFE
675 5020A07	Tube Cutter	1	Plastic with metal blade
675 5020A04	Sample Bottle 180ml	1	PFA
675 5020A98	Tube Weight	1	Stainless Steel

*Sampling Systems offer a PTFE bottle adapter that can be custom made to suit any bottle with a screw thread

Recommended Storage Conditions:

Recommended Spares:

Dry and ambient temperature

5020T-05	PTFE Tubing – 5m length
5020T-100	PTFE Tubing – 100m length
5020T-100D	PTFE Tubing – 100m length & documentation**
5020T-03	Sample Bottle 90ml

** Documentation Package includes a certification of FDA compliance, date of manufacture.

Important Note: If the user continues to operate the pump once the bottle is full the liquid will be drawn into the body of the pump. The pump will then be contaminated and the liquid will leak out.

QAQC LAB 589 Rappahannock Drive WHITE STONE VA 22578 TEL (866) 244-1578