

Part No. 675 8060H510S


Description

Disposable Powder Thief – 500mm – 10ml Tip

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004



Method of Construction:

Ultrasonic welding

Nominal Length:

511mm

Diameter of Outer Tube:

21mm

Nominal Weight of Sampler:

93g

Assembly & Packing Environment:

Class 100,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE bag)

Method of Sterilisation:

Gamma Irradiation (25-45kGy)

Number of Samplers per Box:

20

BSE/TSE Free:

Yes

Recommended Storage Conditions:

Dry and ambient temperature

Note: This is a single-use, disposable sampler. Repeated use will lead to failure.

QAQC LAB WHITE STONE VA 22578 TEL (804) 318-3686

www.qclabequipment.com