Part No. Î Ï Í Á8060H520

QAQE LAB


Description

Disposable Powder Thief - 500mm - 20ml Tip

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004



Method of Construction: Ultrasonic welding

Nominal Length: 511mm

Diameter of Outer Tube: 21mm

Nominal Weight of Sampler: 84g

Assembly & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: NA

Number of Samplers per Box: 20

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

Note: This is a single-use, disposable sampler. Repeated use will lead to failure.

QAQC LAB WHITE STONE VA 22578 TEL (804) 318-3686

www.qclabequipment.com